Series 4560 Mini Reactor Systems

Series Number:

Type: **Mini**

Stand: Bench Top

Vessel Mounting: Moveable or Fixed Head

Sizes, mL: **100 - 600**

Standard Temperature MAWP Rating, psi (bar): **3000 (200)**

High Temperature (HT) MAWP Rating psi (bar): **2000 (138)**

Standard Maximum Operating Temp., °C: 225 w/ FKM 0-ring 275 w/ FFKM 0-ring 350 w/ PTFE Flat Gasket

High Temperature (HT) Maximum Operating Temperature, °C: **500 w/ FG Flat Gasket**

hese are the most popular of all Parr Stirred Reactors. Although they are called "Mini" reactors, they offer a range of sizes large enough to work with significant sample sizes, yet small enough to be handled with ease by all operators. They are made in both fixed head and moveable vessel styles, with a choice of either a self sealing O-ring closure with no cap screws for working temperatures up to 225 °C, or with a flat, PTFE gasket for higher temperatures up to 350 °C sealed with six cap screws in a split-ring cover clamp.

Choosing the high temperature option (HT) boosts the maximum temperature to 500 °C, but de-rates the MAWP to 2000 psi (138 bar). Standard Mini reactors can be converted to high temperature reactors (500 °C max temperature and 2000 psi MAWP) by changing the head openings to cone connections, replacing the valves with high temperature valves, replacing the gasket with a graphoil Flat Gasket, replacing the heater with a ceramic fiber heater, and replacing the split ring. Contact Parr for details.

Although internal and external space is limited in these small vessels, gas entrainment impellers, catalyst baskets, condensers and other options are available. All reactors in this series can also be easily converted from one size to another by simply substituting a longer or shorter cylinder with the corresponding internal fittings and heaters.

The support system for these Mini Reactors is designed specifically to provide stability at stirring speeds up to 1700 rpm, in a compact mounting small enough to fit into a laboratory hood. The support system can also be adapted to accept any of the smaller vessels from the 4590 Micro Reactor Series.



4566 Mini Reactor, 300 mL, Fixed Head, Quick-release Split Ring, and a 4848 Controller shown with optional Expansion Modules.



4566 with Aluminum Block Heater, Bottom Drain Valve, and a 4848 Controller shown with optional Expansion Modules.

Parr

Series 4560 Pressure Reacto	or System	n Specific	ations	-		-				
						vithin series.				
Model Number	4561	4562	4563	4564	4565	4566	4567	4568	4566B	45660
Sizes, mL	300	450	600	160	100	300	450	600	160	100
Maximum Pressure (MAWP)			1		3000 psi	(200 bar)				
HT Max. Pressure (MAWP)		1			2000 psi	(138 bar)				
Maximum Temperature										
with FKM 0-ring					225	5 °C				
with FFKM 0-ring					275	5 °C		-		
with PTFE Flat Gasket				1	350	0°C				
with FG Flat Gasket (HT)		500 °C		N	A		500 °C		N	А
Vessel Style			Moveable					Fixed Head		
Reactor Mounting					Benc	h Top				
Closure					Split-Ring (6	Cap Screws	.)			
Valve Connections					1/8" NP	PT Male				
Magnetic Stirrer, Model No.					A112	OHC6				
Maximum Torque					16 Inch	-Pounds				
Impeller(s), 4-Blade, 1.38" dia.	1	1	2	1	1	1	1	2	1	1
Pressure Gage, Size		1		1	3.5 ir	iches		,		
Range, Standard Pressure			•		0-2000 ps	i (140 bar)				
Range, High Pressure					0-3000 ps	i (200 bar)				
Temperature Measurement	Fixed Thermocouple (Thermowell for special alloys)									
Cooling Coil		Included		Not In	cluded		Included		Not In	cluded
Style					Singl	e Loop				
Bottom Drain Valve	1/4" NPT (Not available on HT Models.)									
Heater Style		Mantle			ip-on		Mantle	-	Clam	ip-on
Heater Power, Watts	510	590	780		00	510	590	780		00
Heater Style w/ HT Option	(Ceramic Fibe	er		/A	(Ceramic Fibe	er		/A
Heater Power, Watts	800	800	1100			800	1000	1200		
Stirrer Motor					1/8 hn vari	able speed				
Electrical Supply					.,					
Volts, AC					115 c	or 230				
Maximum Load, amps, 115 / 230	115 or 230 10 / 5 (14 / 7 High Temperature)									
Vessel Dimensions		1		10,	o (117 7 Hig	jii romporuu	ano,			
Inside Diameter, inches	2.5	2.5	2.5	2.5	2	2.5	2.5	2.5	2.5	2
Inside Depth, inches	4	6	8	2.3	2	4	6	8	2.3	2
Weight of Vessel, pounds	10	13	16	8	9	12	15	18	8	9
Reactor Dimensions	10	10	10	0		12	13	10		
Width, inches w/o Controller					1	2				
Depth, inches						8				
Height, inches*	28	32	36	28	28	28	32	36	28	28
neight, molles										
Woight nounde w/ Controller	60	60	23	0.0	60	62	29	60	0.0	60
Weight, pounds w/ Controller Spare Parts Kit	60	63	66	60	60	62 69M	65	68	60	60

Other options available. See Ordering Guide, visit www.parrinst.com, or call for more information.

1 - 8 0 0 - 8 7 2 - 7 7 2 0

Series 4560 Mini Reactor Systems

INDEX TO OPTIONS

OPTIONS	PAGE
Certification	9
Materials of Construction	10
Magnetic Drive	14
Gaskets & Seals	18
Stirrer Motor	104
Rupture Disc	110
Pressure Gage	111
Cooling Coil	115
Bottom Drain Valve	116
4848 Reactor Control	91
4838 Temperature Control	94
4871 Process Control	95

CUSTOM OPTIONS

Stirrers	106
Heaters	102
External Valves and Fittings	117
Solids Charging Ports	114
Condensers	109
Electrical Glands	122
Catalyst Addition Device	114
Catalyst Baskets	108
Relief Valves	111
Windows	121
Gear Drives	104
Explosion Proof Options	119
Sample Collection Device	116

ACCESSORIES

Liners	115
Spare Parts Kits	122
Pressure Hoses	118
Check Valves	117
Liquid Pipettes	113
Gas Burettes	113



Model 4563 Mini Reactor, 600 mL, Moveable Vessel, PTFE Flat Gasket Seal, with vessel disassembled, and a 4848 Controller shown with optional Expansion Modules.



4560 Reactor Vessels from left to right, 300 mL, 100 mL, 450 mL, 160 mL, and 600 mL.

Series 4560 Ordering Guide

Par

The Order No. for the Base System is: 456_-T-SS-M-115-VS.12-2000-4848

A composite identification number to be used when ordering a 4560 Series Reactor can be developed by combining individual symbols from the separate sections below. For more information on how to use this ordering guide, please see page 27.

A Base Model			
Model No.	Size	Vessel Style	
4561	300 mL	Moveable	
4562	450 mL	Moveable	
4563	600 mL	Moveable	
4564	160 mL*	Moveable	
4565	100 mL *	Moveable	
4566	300 mL	Fixed Head	
4567	450 mL	Fixed Head	
4568	600 mL	Fixed Head	
4566B	160 mL*	Fixed Head	
4566C	100 mL*	Fixed Head	

* High temperature option not available.

B High Ter	nperature Option
-No Symbol	Standard Configuration
-HT	High Temperature / 500 °C
-No Symbol	Standard Configuration

C Gasket	/ Maximum Temperature
-0V	FKM O-ring / 225 °C
-0K	FFKM O-ring / 275 °C
-T	PTFE Flat Compression Gasket / 350 °C
-FG	Flexible Graphite Gasket / 500 °C

D Material	of Construction	
-SS	T316 Stainless Steel	
-M0	Alloy 400*	
-IN	Alloy 600	
-HB	Alloy B-2*	
-HC	Alloy C-276	
-CS	Alloy 20*	
-TI2	Titanium Grade 2*	
-TI4	Titanium Grade 4*	
-NI	Nickel 200*	
-ZR	Zirconium Grade 702 or 705*	
* Not available	* Not available for HT option	

E Magnetic Stirrer Drive		
-M	General Purpose, 16 in-lb	
-FMD1	Footless, 16 in-lb	

F Stirrer Drive Material of Construction		
-No Symbol	Same Material As Vessel	
-MOC Symbol	Indicate Material of Construction	

G Electrica	I Supply	
-115	115 VAC	
-230	230 VAC	

H Motor O	ption
-VS .12	Variable Speed, 1/8 hp
-VS .25	Variable Speed, 1/4 hp
-XP .25	Explosion Proof Variable Speed, 1/4 hp
-AM .25	Air Motor, 1/4 hp

I Pressure Gage	
-3000	3000 psi / 200 bar
-2000	2000 psi / 137 bar
-1000	1000 psi / 70 bar
-600	600 psi / 40 bar
-200	200 psi / 14 bar
-100	100 psi / 7 bar

J Bottom Drain Valve	
-No Symbol	No Bottom Drain Valve
-BDV	Bottom Drain Valve Installed (Not Available on High Temperature Vessels)

K Controller		
-4848	PID Control, Ramp & Soak Programming, Motor Speed Control, and Data logging with Software	
-A1925E4	Optional RS485 to USB Cable for 4848 Controller	
-4871	Process Controller	

L 4848 Expansion Modules	
-TDM	Tachometer Display Module
-MCM	Tachometer w/Motor Control Module
-PDM	Pressure Display Module
-HTM	High Temperature Cut Off Module
-ETLM	External Temperature Limit Module
-MTM	Motor Torque Module

M 4848 Cooling Control: Solenoid Valve Module		
-SVM1		115 VAC
-SVM2		230 VAC

N Custom Options (List All Desired)	
-AS	Anchor Stirrer
-SB*	Static Catalyst Basket
-CAD*	Catalyst Addition Device
-BF*	Baffles
-GE	Gas Entrainment Stirrer
-RC	Reflux Condenser
-RTC	Reflux Take-off Condenser
-SCP	Solids Charging Port
-ABH*	Aluminum Block Heater
-WJ	Welded Jacket
* 300 mL vessels and larger only	

0 Certifications	
No Certification	
ASME Certification	
European Community Certification/Pressure	
China Certification	
Parr Certification	

P Spare Parts Kit	
-4569M	Spare Parts Kit for 4560 Series