



## **NIRMasteR™** Meat and Meat Products

Production of meat products is still based on a high degree of trust. To ensure high quality and legislation fulfillments it is very important to have access to accurate and meaningful quality parameters. Understanding customer needs in the meat industry, BUCHI has designed the first standalone FT-NIR Spectrometer specifically for that industry – the NIRMasteR™.



## **NIRMaster™** Bring Light into Your Production

Get more analysis information, higher productivity, higher quality and better understanding in meat production with polarization Fourier-Transform Near-Infrared Spectroscopy (FT-NIRS). Rapidly control the quality and nutritional facts of your meat products using BUCHI NIRsolutions. The NIRMaster™ is the first dedicated FT-NIR standalone spectrometer. Its easy-clean design and high ingress protection (IP) meet the specific needs for at-line analysis in meat production.

### Where BUCHI NIRsolutions can help you:

BUCHI NIRsolutions provides dedicated at-line and off-line NIR equipment. These instruments set up at different locations in the production plant provide significant information for meat processing and quality control to improve your business.



### A contribution to your profit

*“Polarization FT-NIR enables for measuring the fat quality of carcasses already in the slaughter plant. In Switzerland fat quality is part of the payment system. In contrast to the conventional iodine value FT-NIR gives additional information about the fatty acid profile. The measurement is done easier, faster and without the use of toxic solvents. Due to the transferable calibrations, labs at various plants do measure the exactly same trait, which is tremendously important for an economically relevant quality criterion.”*  
SUISAG

---

Various parameters of meat products such as moisture, fat, protein, ash, drip loss, collagen, fat score, iodine value and many more can be qualified and quantified. Please refer to our application overview.

---

## 1. Quality control of raw meat and ingredients of meat products

In the slaughterhouse or at the intake:

- Identify and qualify incoming goods (e.g. origin)
- Rapid determination of quality parameters
- Rapidly quantify raw material and components
- Establish a firewall against imposture



## 2. Determination of nutritional aspects and production control

At production site:

- Monitor manufacturing processes
- Optimize your formulation
- Rapidly quantify constituents
- Control component levels (e.g. blending with a specific level of fat)



## 3. Quality control of finished meat products

At packaging area:

- Rapidly release products
- Check production specification labelling
- Remove bottlenecks in production and delivery





## **NIRMaster™ Makes Your Job Easy** The Dedicated Solution for Your Needs

The competitiveness of the meat processing industry forces producers to take strategic actions for quality control and product optimization. BUCHI FT-NIR polarization technology makes this job easier!



### **1. Easy-clean design increases food safety day by day**

The hygienic easy-clean design and high ingress protection (IP) support your routine work in harsh environments. The NIRMaster™ is available in food-grade PMMA housing or stainless steel housing.

### **2. Polarization FT-NIR reduces costs**

The modern and robust NIRMaster™ based on fast and precise polarization FT-NIR technology, provides accurate and reliable results for multiple parameters simultaneously. Improve your production process with this fast and durable system.

### **3. Faster results speed up your decisions**

The NIRMaster™ with its built in PC, preinstalled software and quick-start calibrations can be quickly set up and put into operation. With FT-NIR you can analyze your samples within seconds.

### **4. Reliable hardware for reliable business**

The NIRMaster™ is supported by a twin-lamp safeguard for uninterrupted operation even in cases of lamp failure. Based on robust and accurate polarization technology it is more shock resistant than any other FT-NIR.

### **5. Easy-operation for a lean workflow**

The NIRWare Operator is the easy to use and ergonomic software interface which can be tailored to your specific needs. Inputs are collected easily and the analysis results are displayed in an optimal way. With an optional touch screen and robust accessories, like a barcode scanner or the easy-spin drive, efficiency in daily work can be improved even more.



## 6. Quick-start calibrations to gain time

In cooperation with many of BUCHI's long-term customers we have concentrated data and experience in our pre-calibrations. These powerful quick-start calibrations enable you to efficiently analyze a broad range of your sample constituents.

## 7. Easy calibration transfer preserves competence

BUCHI NIR spectrometers are identical by design. Spectra and calibrations can be exchanged between laboratory and production units - without any calibration standardization. Furthermore, your existing calibrations, even data from other suppliers, can be transferred easily to your new BUCHI FT-NIR spectrometer.

## 8. NIRCAl creates better performance

Use the powerful NIRCAl software to evaluate your production data. Keep an eye on what is going on at your sites. For evaluating or updating a calibration NIRCAl is the ideal tool. The patented Calibration Wizard and the comprehensive NIRCAl Toolbox support even unexperienced users building up powerful calibrations.

## 9. NIR Networking for local and global control

For local data handling the LIMS interface connects your BUCHI NIR to laboratory information or enterprise resource planning software using standard interchange formats. For global players NIRAnywhere is the comprehensive network management tool. Headquarters can automatically backup data and update applications by remote connection. You gain full access and transparency to all your production and QC relevant data.

## 10. NIRMaste<sup>TM</sup> as part of all-embracing NIRSolutions

BUCHI NIRSolutions not only offers leading hardware technology, but also the most comprehensive standard software, applications and reference methods as well as quality service and support to help you work effectively and efficiently. The BUCHI NIRMaste<sup>TM</sup> is an integral part of BUCHI NIRSolutions.



## **NIR Solutions®** Listening, Understanding and Responding

BUCHI is the world leading supplier of dedicated solutions for the meat industry. We provide analytical equipment fitting the needs of our customers. The long history of BUCHI in this special market segment gives our customers the possibility to choose the best solution for their needs. We provide the classical reference methods and NIR for quality control of meat products.



### **Enter the NIR Solutions community**

The BUCHI NIR Solutions community is built by our international customers from the meat processing industry and our global expert network. Based on experience and know-how we have designed and created the BUCHI NIRMaste<sup>TM</sup> together.

### **Global expert network**

BUCHI has developed a global unique network of NIR specialists in meat process control, present in more than 60 countries. Our network is based on shared experiences, fast and flexible teams, intensive trainings and access to all relevant resources. These synergies are key features for our unique support structure.

### **Global service network & training**

Our excellence is global and we have local services to fulfill your requirements. We offer service contracts, preventive maintenance protocols, IQ/OQ, on-site service and training. We guarantee fast response, phone and web-based support. For our key customers we offer tailored solutions, services and support.

### **Reference methods & feasibility studies**

BUCHI has gained high expertise and market share with relevant reference methods like Kjeldahl or extraction, which are official methods in meat process control. The accuracy and precision of our NIR applications is based on the knowledge of both, NIR and reference methods. We support you with feasibility studies and professional competence.




---

## Advantages of BUCHI NIR Spectrometers:

Compared to classical analytical methods the NIRMaste<sup>TM</sup>

- Enables you to have faster analytical results
  - Speeds up decision making
  - Grants you low analysis costs per sample
  - Measures several parameters simultaneously
  - Improves quality and optimizes overall profit
- 

---

## Advantages of polarization FT-NIR over any other NIR technology

- Most accurate and reproducible NIR technology
  - The polarization interferometer is insensitive to vibration
  - Increase of information content and application capabilities
  - Calibration transfer without any standardization
  - Twin-lamp safeguard for uninterrupted operation
- 

## BUCHI NIR Solutions covers all your needs - from Laboratory to Production Environment

Using the same technology as The NIRMaste<sup>TM</sup> the NIRFlex N-500 is the flexible solution for your laboratory. Swap between different measuring cells – fitting your demands. Our customized solutions will make your BUCHI FT-NIR spectrometer become more versatile.



# NIRMaster™ The Complete Solution



## Performance

- Measurement of several parameters simultaneously
- Low analysis costs per sample



## Durability

- Robust technology
- High ingress protection (IP)
- Standardization-free calibrations



## Speed

- Quick-start calibrations
- Faster analytical results
- Speed up decision making



## Reliability

- Accurate and robust polarization FT-NIR
- Twin-lamp safeguard



## Simplicity

- Built-in PC
- Touch screen use
- Easy-clean design



## Competence

- Comprehensive applications
- Global expert network
- Best service and support

## Quality in your hands

### Find your BUCHI Contacts:

**BÜCHI Labortechnik AG**  
CH – 9230 Flawil 1  
T +41 71 394 63 63  
F +41 71 394 65 65  
buchi@buchi.com  
www.buchi.com

**BUCHI UK Ltd.**  
GB – Oldham OL9 9QL  
T +44 161 633 1000  
F +44 161 633 1007  
uk@buchi.com  
www.buchi.co.uk

**BUCHI Hong Kong Ltd.**  
HK – Central  
T +852 2389 2772  
F +852 2389 2774  
china@buchi.com  
www.buchi.com.cn

**Nihon BUCHI K.K.**  
JP – Tokyo 110-0008  
T +81 3 3821 4777  
F +81 3 3821 4555  
nihon@buchi.com  
www.nihon-buchi.jp

**BUCHI Korea Inc**  
KR – Seoul 153-782  
T +82 2 6718 7500  
F +82 2 6718 7599  
korea@buchi.com  
www.buchi.com

**BÜCHI Labortechnik GmbH**  
DE – 45127 Essen  
Freecall 0800 414 0 414  
T +49 201 747 490  
F +49 201 237 082  
deutschland@buchi.com  
www.buechigmbh.de

**BÜCHI Labortechnik GmbH**  
NL – 3340 AC  
Hendrik-Ido-Ambacht  
T +31 78 684 94 29  
F +31 78 684 94 30  
netherlands@buchi.com  
www.buchi.nl

**BUCHI Shanghai**  
CN – 200052 Shanghai  
T +86 21 6280 3366  
F +86 21 5230 8821  
china@buchi.com  
www.buchi.com.cn

**BUCHI India Private Ltd.**  
IN – Mumbai 400 055  
T +91 22 667 18983 / 84 / 85  
F +91 22 667 18986  
india@buchi.com  
www.buchi.in

**BUCHI Canada Ltd.**  
CA – Thornhill, Ontario L4J 6Z2  
T +1 416 277 7407  
F +1 905 764 5139  
canada@buchi.com  
www.mybuchi.com

**BUCHI Sarl**  
FR – 94656 Rungis Cedex  
T +33 1 56 70 62 50  
F +33 1 46 86 00 31  
france@buchi.com  
www.buchi.fr

**BUCHI Italia s.r.l.**  
IT – 20090 Assago (MI)  
T +39 02 824 50 11  
F +39 02 57 51 28 55  
italia@buchi.com  
www.buchi.it

**BUCHI (Thailand) Ltd.**  
TH – Bangkok 10600  
T +66 2 862 08 51  
F +66 2 862 08 54  
bacc@buchi.com  
www.buchi.com

**PT. BUCHI Indonesia**  
ID – Tangerang 15321  
T +62 21 537 62 16  
F +62 21 537 62 17  
indonesia@buchi.com  
www.buchi.com

**BUCHI Corporation**  
US – New Castle,  
Delaware 19720  
Toll Free: +1 877 692 8244  
T +1 302 652 3000  
F +1 302 652 8777  
us-sales@buchi.com  
www.mybuchi.com

We are represented by more than 100 distribution partners worldwide.  
Find your local representative at: [www.buchi.com](http://www.buchi.com)

