

Non-Stirred Pressure Vessels

General Purpose Non-stirred Reactors and Pressure Vessels

Parr offers a broad selection of non-stirred pressure vessels in convenient styles, sizes and pressure ranges for many laboratory uses. Recent additions to this line have included new sizes, new temperature and pressure combinations and new self-sealing closures, all designed to provide convenient vessels for laboratory procedures which must be performed under pressure at elevated temperatures. For more information, please read our 'Selection Guide to Parr Pressure Vessels.'



Series 4600 - 4620 GP 4.00 Inch Diameter

Sizes: 1000mL, 2000mL



Series 4632 Pressure Extraction Vessel 4.00 Inch Diameter

Size: 2000 mL



Series 4650 2.50 Inch Inside Diameter HP/HT

Sizes: 250mL, 500mL, 1000mL



Series 4660 6.00 Inch Inside Diameter GP

Sizes: 3.75L, 7.5L



Series 4670 5.50 Inch Inside Diameter HP/HT

Sizes: 3750mL, 5800mL



Series 4680 3.75 Inch Inside Diameter HP/HT

Sizes: 1000mL, 1800mL



Series 4676 9.50 Inch Diameter GP

Sizes: 10000, 20000mL



Series 4750 1.50 Inch Inside diameter GP

Sizes: 125mL, 200mL



Series 4760 - 4774 2.50 Inch Inside Diameter GP

Sizes: 160mL, 600mL



Series 4790 Micro GP

Sizes: 25 mL, 50mL, 100 mL



Series 4703 - 4714 Screw Cap Vessels 1.00 Inch Inside Diameter

Sizes: 22mL, 45mL



Series 4740 1.00 Inch Inside Diameter HP/HT

Sizes: 21mL, 71mL