

SpeedExtractor E-914 / E-916 Combine maximum speed and throughput

The SpeedExtractor is your best solution for fast Pressurized Solvent Extraction (PSE). Increase productivity by processing up to 6 samples in parallel. Strea^{mL}ine the workflow of your sample preparation thanks to ease of sample loading and ready-to-use extract collection.



Unrivaled throughput Extract 6 samples in 30 min and free up your time for other tasks

Ease of operation

Automated sealing reduces operator intervention and increases process safety





Low running costs

Combine extract collection with evaporation, minimize solvent use and avoid costly consumable seals



Unrivaled throughput

- · 6 times faster than other PSE instruments thanks to parallel extraction concept
- Extract up to 96 samples in an 8 hour shift
- · Accelerate the extraction process and facilitate final results within one working day
- \cdot Reduce number of replicates due to identical extraction conditions for up to 6 samples



Ease of operation

- \cdot Unique design of extraction cells allows for ease of sample loading and cell assembly
- \cdot High level of automation saves time-consuming preparation and increases safety
- \cdot Quick, fail-safe and reproducible operation due to assembly of only 4 sample cell parts
- \cdot The automated sealing eliminates sources of error and guarantees safety standard



Low running costs

- Complementary workflow and synergy with parallel evaporation and concentration
 Collection vessels with appendix make extract transfer obsolete
- · Reduced solvent consumption and minimal need for consumables give lower running costs
- \cdot Save energy and costs with ECO mode for automated heater control

Complete your portfolio



Syncore® Analyst Maximize multiple sample workup efficency



Multivapor P-6 / P-12 Efficient evaporation for multiple samples



Sepacore® Easy Extract purification system Clean up of complex extracts



Hydrolysis Unit E-416 / B-411 Safe and fast acid digestion